

For Research Use Only

TMEM150 Recombinant monoclonal antibody

Catalog Number: 87578-1-RR



Basic Information

| | | |
|--|--|---|
| Catalog Number: 87578-1-RR | GenBank Accession Number: BC050466 | Purification Method: Protein A purification |
| Size: 100ul , Concentration: 1000 µg/ml by Nanodrop; | GeneID (NCBI): 129303 | CloneNo.: 252959H6 |
| Source: Rabbit | UNIPROT ID: Q86TG1 | Recommended Dilutions: WB: 1:500-1:2000 |
| Isotype: IgG | Full Name: transmembrane protein 150 | |
| Immunogen Catalog Number: AG36287 | Calculated MW: 271 aa, 29 kDa | |
| | Observed MW: 29 kDa | |

Applications

| | |
|--|---|
| Tested Applications: WB, ELISA | Positive Controls: WB : Jurkat cells, Neuro-2a cells, PC-12 cells |
| Species Specificity: human, mouse, rat | |

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87578-1-RR (TMEM150 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.