

For Research Use Only

- PEG Recombinant monoclonal antibody, PBS Only

Catalog Number:87604-1-PBS



Basic Information

Catalog Number:

87604-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

240925A1

Applications

Tested Applications:

ELISA, Indirect ELISA

Species Specificity:

-

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



96-well plates were coated with Peptid PEG12, mPEG20K-Biotin, Biotin-PEG2K-Biotin. The plate was blocked with 1% BSA blocking buffer, then washed, and incubated with 100 μ L per well of PEG antibody 87604-1-RR and serially diluted with a starting concentration of 2 μ g/mL (1 hr, 37°C). The plate was washed and incubated with 100 μ L per well of an HRP-conjugated goat anti-Rabbit IgG secondary antibody at a dilution of 1:10,000 (1 hr, 37°C). Detection was

