

For Research Use Only

Mouse Abcg2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87614-2-PBS



Basic Information

Catalog Number: 87614-2-PBS	GenBank Accession Number: BC053730	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 26357	CloneNo.: 252770G12
Source: Rabbit	UNIPROT ID: Q7TMS5	
Isotype: IgG	Full Name: ATP-binding cassette, sub-family G (WHITE), member 2	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Product Information

87614-2-PBS targets Abcg2 as part of a matched antibody pair:

MP03130-1: 87614-3-PBS capture and 87614-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BCRP/ABCG2 concentration was determined to be 8.7 ng/mL in mouse liver tissue extract based on a 8.4 mg/mL extract load, mouse eye tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03130-1, MOUSE BCRP/ABCG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87614-3-PBS. Detection antibody: 87614-2-PBS. Standard: SY02955. Range: 0.156-20 ng/mL