

For Research Use Only

# TEKT2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87634-2-PBS



## Basic Information

<b>Catalog Number:</b> 87634-2-PBS	<b>GenBank Accession Number:</b> BC035620	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 27285	<b>CloneNo.:</b> 252977G12
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9UIF3	
<b>Isotype:</b> IgG	<b>Full Name:</b> tektin 2 (testicular)	
<b>Immunogen Catalog Number:</b> AG4457	<b>Calculated MW:</b> 430 aa, 50 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

87634-2-PBS targets TEKT2 as part of a matched antibody pair:

MP03140-1: 87634-3-PBS capture and 87634-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean TEKT2 concentration was determined to be 49.2 ng/mL in mouse testis tissue extract based on a 4.2 mg/mL extract load, 2.1 ng/mL in mouse lung tissue extract based on a 5.0 mg/mL extract load.



Cytometric bead array standard curve of MP03140-1, TEKT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87634-3-PBS. Detection antibody: 87634-2-PBS. Standard: Ag4457. Range: 0.625-80 ng/mL.