

For Research Use Only

USF2 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87636-1-PBS



Basic Information

Catalog Number:

87636-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9866

GenBank Accession Number:

BC049821

GeneID (NCBI):

7392

UNIPROT ID:

Q15853

Full Name:

upstream transcription factor 2, c-fos interacting

Calculated MW:

346 aa, 37 kDa

Observed MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

252980B1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87636-1-RR (USF2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87636-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87636-1-RR against Human USF2 were performed. The affinity constant is below 1 pM.