

For Research Use Only

NCAPG2 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87645-1-PBS



Basic Information

Catalog Number: 87645-1-PBS	GenBank Accession Number: BC043404	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 54892	CloneNo.: 253099H6
Source: Rabbit	UNIPROT ID: Q86XI2	
Isotype: IgG	Full Name: non-SMC condensin II complex, subunit G2	
Immunogen Catalog Number: AG34991	Calculated MW: 1143 aa, 131 kDa	
	Observed MW: 130 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

NCAPG2 (non-SMC condensin II complex subunit G2) encodes a protein belonging to the non-SMC (structural maintenance protein) subunit family of the Condensin II complex that plays a key role in chromosome condensation and segregation during mitosis. NCAPG2 is positively correlated with cancer stem cell (CSC) markers, and aberrant NCAPG2 expression contributes to tumorigenesis in several human malignancies, such as prostate cancer and low-grade gliomas.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87645-1-RR (NCAPG2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87645-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using NCAPG2 antibody (87645-1-RR, Clone: 253099H6) at dilution of 1:500 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 87645-1-PBS in a different storage buffer formulation.