

For Research Use Only

ZNF207 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87649-1-PBS



Basic Information

Catalog Number: 87649-1-PBS	GenBank Accession Number: BC002372	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7756	CloneNo.: 253049C10
Source: Rabbit	UNIPROT ID: O43670	
Isotype: IgG	Full Name: zinc finger protein 207	
Immunogen Catalog Number: AG40360	Calculated MW: 51 kDa	
	Observed MW: 50-57 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87649-1-RR (ZNF207 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87649-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87649-1-RR against Human ZNF207 were performed. The affinity constant is 0.478 nM.