

For Research Use Only

ZNF207 Recombinant monoclonal antibody

Catalog Number: 87649-1-RR



Basic Information

Catalog Number:

87649-1-RR

Size:

100ul, Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG40360

GenBank Accession Number:

BC002372

GeneID (NCBI):

7756

UNIPROT ID:

O43670

Full Name:

zinc finger protein 207

Calculated MW:

51 kDa

Observed MW:

50-57 kDa

Purification Method:

Protein A purification

CloneNo.:

253049C10

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HuH-7 cells, A549 cells, HepG2 cells, Neuro-2a cells, NIH/3T3 cells, 3T3-L1 cells, C6 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

