

For Research Use Only

PAX9 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87678-2-PBS



Basic Information

Catalog Number: 87678-2-PBS	GenBank Accession Number: BC001159	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5083	CloneNo.: 253236E9
Source: Rabbit	UNIPROT ID: P55771	
Isotype: IgG	Full Name: paired box 9	
Immunogen Catalog Number: AG25381	Calculated MW: 341 aa, 36 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

87678-2-PBS targets PAX9 as part of a matched antibody pair:

MP03150-1: 87678-2-PBS capture and 87678-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PAX9 concentration was determined to be 9.6 ng/mL in A2780 cell extract based on a 1.2 mg/mL extract load, SH-SY5Y cell extract were served as a negative control.



Cytometric bead array standard curve of MP03150-1, PAX9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87678-2-PBS. Detection antibody: 87678-1-PBS. Standard: Ag25381. Range: 0.156-20 ng/mL.