

For Research Use Only

CD27/TNFRSF7 Recombinant monoclonal antibody

Catalog Number:87691-3-RR



Basic Information

Catalog Number: 87691-3-RR	GenBank Accession Number: NM_001024335.1	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 500318	CloneNo.: 253071D4
Source: Rabbit	UNIPROT ID: Q501W2	Recommended Dilutions: WB: 1:2000-1:10000
Isotype: IgG	Full Name: CD27 molecule	
Immunogen Catalog Number: EG5288	Calculated MW: 28 kDa	
	Observed MW: 50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : rat spleen tissue, rat kidney tissue, rat liver tissue
Species Specificity: rat	

Background Information

Rat CD27 (TNFRSF7) is a member of the tumor necrosis factor receptor superfamily (TNFRSF). It is a type I transmembrane glycoprotein primarily expressed on the surface of T lymphocytes, B cells, and natural killer (NK) cells. As an important T cell co-stimulatory molecule, CD27 binds to its ligand, CD70. This interaction recruits TRAF2/5, activating the NF-κB and JNK signaling pathways, which promotes T cell proliferation, differentiation, and survival while inhibiting apoptosis. The expression level of CD27 serves as a crucial marker for T cell activation and memory status, and soluble CD27 (sCD27) can act as a clinical marker of immune activation. In cancer immunotherapy, CD27 is a significant target. Agonistic antibodies, such as Varlilumab, can enhance anti-tumor T cell responses and have entered clinical trials.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87691-3-RR (CD27/TNFRSF7 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.