

For Research Use Only

MARCS Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87693-3-PBS



Basic Information

Catalog Number: 87693-3-PBS	GenBank Accession Number: NM_002356.7	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4082	CloneNo.: 253128D11
Source: Rabbit	UNIPROT ID: P29966	
Isotype: IgG	Full Name: myristoylated alanine-rich protein kinase C substrate	
Immunogen Catalog Number: EG6453	Calculated MW: 32 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

87693-3-PBS targets MARCS as part of a matched antibody pair:

MP03172-1: 87693-3-PBS capture and 87693-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MARCKS concentration was determined to be 83.8 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 32.3 ng/mL in SH-SY5Y cell extract based on a 1.5 mg/mL extract load, 8.3 ng/mL in A-204 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03172-1, MARCKS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87693-3-PBS. Detection antibody: 87693-2-PBS. Standard: Eg6453. Range: 0.312-40 ng/mL.