

For Research Use Only

TSPAN8 Recombinant monoclonal antibody

Catalog Number:87700-3-RR



Basic Information

Catalog Number: 87700-3-RR	GenBank Accession Number: NM_004616.2	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 7103	CloneNo.: 252991H1
Source: Rabbit	UNIPROT ID: P19075	Recommended Dilutions: WB: 1:5000-1:50000 IHC: 1:250-1:1000
Isotype: IgG	Full Name: tetraspanin 8	
Immunogen Catalog Number: EG7574	Calculated MW: 26 kDa	
	Observed MW: 27 kDa	

Applications

Tested Applications:
WB, IHC, ELISA

Species Specificity:
human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : Calu-3 cells, Ramos cells, LoVo cells

IHC : human stomach cancer tissue,

Background Information

TSPAN8 is a membrane protein of the tetraspanin superfamily, which are characterized by the presence of four conserved transmembrane regions. Tetraspanins have been involved in diverse processes such as cell activation and proliferation, adhesion and motility, differentiation, and cancer. TSPAN8 and CD151 coordinately promote metastasis, where Tspan8 overrides the adhesive features of CD151 by recruiting integrins out of adhesion into motility promoting complexes (PMID: 23683890). TSPAN8 contributes to molecular pathways of exosome-induced endothelial cell activation (PMID: 20124479).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87700-3-RR (TSPAN8 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 87700-3-RR (TSPAN8 antibody) at dilution of 1:500 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).