

For Research Use Only

CCDC80 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87731-2-PBS



Basic Information

Catalog Number: 87731-2-PBS	GenBank Accession Number: BC116178	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 151887	CloneNo.: 253134D6
Source: Rabbit	UNIPROT ID: Q76M96	
Isotype: IgG	Full Name: coiled-coil domain containing 80	
Immunogen Catalog Number: EG6754	Calculated MW: 950 aa, 108 kDa	
	Observed MW: 110-180 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Coiled-coil domain-containing protein 80 (CCDC80) is a member of the coiled-coil domain-containing protein family, which is mainly expressed in vascular smooth muscle cells. It plays an integral role in the regulation of cardiovascular remodeling and homeostasis (PMID: 23831687; 30439364; 27375481). In SDS PAGE, CCDC80 can be detected as double bands (PMID: 40278823; 35843176).

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Caki-1 cells were subjected to SDS PAGE followed by western blot with 87731-2-RR (CCDC80 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87731-2-PBS in a different storage buffer formulation.