

For Research Use Only

FRS3 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87810-1-PBS



Basic Information

Catalog Number:

87810-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2896

GenBank Accession Number:

BC010611

GeneID (NCBI):

10817

UNIPROT ID:

O43559

Full Name:

fibroblast growth factor receptor substrate 3

Calculated MW:

492 aa, 54 kDa

Observed MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

253396B5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com


This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87810-1-RR (FRS3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87810-1-PBS in a different storage buffer formulation.

Various lysates were subjected to SDS PAGE followed by western blot with 87810-1-RR (FRS3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87810-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87810-1-RR against Human FRS3 were performed. The affinity constant is 0.955 nM.