

For Research Use Only

# FRS3 Recombinant monoclonal antibody

Catalog Number: 87810-1-RR



## Basic Information

Catalog Number:

87810-1-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2896

GenBank Accession Number:

BC010611

GeneID (NCBI):

10817

UNIPROT ID:

O43559

Full Name:

fibroblast growth factor receptor substrate 3

Calculated MW:

492 aa, 54 kDa

Observed MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

253396B5

Recommended Dilutions:

WB: 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : DU 145 cells, mouse brain tissue, PC-3 cells, rat brain tissue

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87810-1-RR (FRS3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.

Various lysates were subjected to SDS PAGE followed by western blot with 87810-1-RR (FRS3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.