

For Research Use Only

IFNA8 Recombinant monoclonal antibody

Catalog Number: 87827-1-RR



Basic Information

Catalog Number: 87827-1-RR	GenBank Accession Number: BC104830	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 3445	CloneNo.: 252963A7
Source: Rabbit	UNIPROT ID: P32881	Recommended Dilutions: WB: 1:2000-1:10000
Isotype: IgG	Full Name: interferon, alpha 8	
Immunogen Catalog Number: EG4954	Observed MW: 25 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Transfected with IFNA8 HEK-293T cells,
Species Specificity: human	

Background Information

IFNA8 (Interferon alpha-8) is a member of the type I interferon family, primarily produced by macrophages and virus-infected cells, and possesses broad-spectrum antiviral activity. This protein exerts its effects by binding to the type I interferon receptor, activating the JAK-STAT signaling pathway, and inducing the expression of various interferon-stimulated genes (ISGs). Studies have shown that the antiviral activity of IFNA8 is significantly higher than that of other subtypes such as IFNA1, playing a crucial role in defending against infections like vesicular stomatitis virus. Furthermore, the transcription of IFNA8 is regulated by the transcription factor Oct-1, and single nucleotide polymorphisms in its promoter region are associated with the prognosis of glioma patients.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87827-1-RR (IFNA8 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 87827-1-RR against Human IFNA8 were performed. The affinity constant is 0.54 nM.