

For Research Use Only

FCRL2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87834-3-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 87834-3-PBS | GenBank Accession Number: NML_030764.3 | Purification Method: Protein A purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 79368 | CloneNo.: 253006C2 |
| Source: Rabbit | UNIPROT ID: Q96LA5-1 | |
| Isotype: IgG | Full Name: Fc receptor-like 2 | |
| Immunogen Catalog Number: EG4773 | Calculated MW: 56 kDa | |

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

87834-3-PBS targets FCRL2 as part of a matched antibody pair:

MP03222-2: 87834-1-PBS capture and 87834-3-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean FCRL2 concentration was determined to be 6.75 ng/mL in Daudi cell extract based on a 1.2 mg/mL extract load and 5.41 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load.



Plasma of six individual healthy human donors was measured. The human FCRL2 concentration of detected samples was determined to be 1,294.6 pg/mL with a range of 980.1-1,641.9 pg/mL.



Sandwich ELISA standard curve of MP03222-2, Human FCRL2 Recombinant Matched Antibody Pair - PBS only. 87834-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg4773. 87834-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL