

For Research Use Only

TNFRSF6B Recombinant monoclonal antibody

Catalog Number: 87840-1-RR



Basic Information

Catalog Number: 87840-1-RR	GenBank Accession Number: NM_003823.3	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 8771	CloneNo.: 253457A6
Source: Rabbit	UNIPROT ID: O95407	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: tumor necrosis factor receptor superfamily, member 6b, decoy	
Immunogen Catalog Number: EG3110	Calculated MW: 33 kDa	
	Observed MW: 35 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls:
Species Specificity: human	WB : ACHN cells, SW480 cells, HCT 116 cells, SW 1990 cells

Background Information

The tumor necrosis factor receptor super-family member 6b (TNFRSF6B) is over-expressed in various human cancers. TNFRSF6B mRNA and protein expression in the tumor tissues negatively correlated with apoptosis. Cell proliferation was decreased, cell cycle was arrested in G1/S-phase, apoptosis was increased, and migration ability was inhibited by anti-TNFRSF6B McAb in vitro. Anti-TNFRSF6B McAb could be useful to suppress proliferation and induce apoptosis in HCC. Thus, TNFRSF6B might be a critical, targeted therapy strategy for HCC (PMID: 20591579).

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87840-1-RR (TNFRSF6B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



ACHN cells were subjected to SDS PAGE followed by western blot with 87840-1-RR (TNFRSF6B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 87840-1-RR against Mouse TNFRSF6B were performed. The affinity constant is 0.131 nM.