

For Research Use Only

CDR2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87868-1-PBS



Basic Information

Catalog Number: 87868-1-PBS	GenBank Accession Number: BC017503	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 1039	CloneNo.: 253563C10
Source: Rabbit	UNIPROT ID: Q01850	
Isotype: IgG	Full Name: cerebellar degeneration-related protein 2, 62kDa	
Immunogen Catalog Number: AG24100	Calculated MW: 454 aa, 52 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

87868-1-PBS targets CDR2 as part of a matched antibody pair:

MP03234-1: 87868-1-PBS capture and 87868-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CDR2 concentration was determined to be 374.0 ng/mL in HepG2 cell extract based on a 1.5 mg/mL extract load, 100.6 ng/mL in A2780 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03234-1, CDR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87868-1-PBS. Detection antibody: 87868-2-PBS. Standard: Ag24100. Range: 0.781-100 ng/mL