

For Research Use Only

Phospho-ALK (Tyr1278) Recombinant monoclonal antibody

Catalog Number: 87874-1-RR



Basic Information

Catalog Number:

87874-1-RR

Size:

100ul, Concentration: 1000 µg/ml by 238

Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_004304

GeneID (NCBI):

238

UNIPROT ID:

Q9UM73

Full Name:

anaplastic lymphoma receptor

tyrosine kinase

Calculated MW:

176 kDa

Observed MW:

80 kDa, 300 kDa

Purification Method:

Protein A purification

CloneNo.:

253445D5

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : pervanadate treated Karpas-299 cells,

Background Information

Anaplastic lymphoma kinase was first described in 1994 as the NPM-ALK fusion protein that is expressed in the majority of anaplastic large-cell lymphomas. About 2/3rds of anaplastic large-cell lymphomas possess a balanced chromosomal translocation in which the entire nucleophosmin (NPM) gene on chromosome 5 is fused to the 3' portion of the ALK gene (including the entire intracellular segment with its protein kinase domain) on chromosome 2. This oncogenic ALK protein kinase is a chimeric protein created by a translocation between chromosomes (2;5) (p23;q35) that generates the NPM-ALK fusion protein. Phospho-ALK (Tyr1278) is a hallmark of ALK-driven malignancies. (PMID: 30400214)

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Non-treated and pervanadate treated Karpas-299 cells were subjected to SDS PAGE followed by western blot with 87874-1-RR (Phospho-ALK (Tyr1278) antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. The membrane was stripped and re-blotted with Vinculin (86930-1-RR) antibody as loading control.