

For Research Use Only

NFATC2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87895-3-PBS



Basic Information

Catalog Number: 87895-3-PBS	GenBank Accession Number: BC136418	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4773	CloneNo.: 253473B2
Source: Rabbit	UNIPROT ID: Q13469	
Isotype: IgG	Full Name: nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 2	
Immunogen Catalog Number: AG16849	Calculated MW: 925 aa, 100 kDa	

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

87895-3-PBS targets NFATC2 as part of a matched antibody pair:

MP03256-1: 87895-3-PBS capture and 87895-2-PBS detection (validated in Cytometric bead array)

MP03256-2: 87895-3-PBS capture and 87895-4-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




The mean NFATC2 concentration was determined to be 5,167.3 pg/mL in MDA-MB-453 cell extract based on a 1.2 mg/mL extract load and 2,012.9 pg/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP03256-2, Human NFATC2 Recombinant Matched Antibody Pair - PBS only. 87895-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag16849. 87895-4-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean NFATC2 concentration was determined to be 2.5 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 27.7 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, RAW267.4 cell extract were served as a negative control.



Cytometric bead array standard curve of MP03256-1, NFATC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87895-3-PBS. Detection antibody: 87895-2-PBS. Standard: Ag16849. Range: 0.156-20 ng/mL.