

For Research Use Only

EPHA2 Recombinant monoclonal antibody

Catalog Number: 87904-1-RR



Basic Information

Catalog Number: 87904-1-RR	GenBank Accession Number: NM_010139.3	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 13836	CloneNo.: 253572H2
Source: Rabbit	UNIPROT ID: Q03145	Recommended Dilutions: WB: 1:2000-1:10000 IF/ICC: 1:200-1:800
Isotype: IgG	Full Name: Eph receptor A2	
Immunogen Catalog Number: EG8778	Calculated MW: 109 kDa	
	Observed MW: 130 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls:
Species Specificity: mouse	WB : L-929 cells, NIH/3T3 cells IF/ICC : NIH/3T3 cells,

Background Information

EPH receptor A2 (EPHA2) is a receptor tyrosine kinase (RTK) of the Eph receptor family that primarily binds glycosylphosphatidylinositol (GPI)-anchored ephrin-A ligands, especially ephrin-A1 (EFNA1), to mediate contact-dependent cell-cell signaling. In mice, Epha2 is broadly expressed in epithelial, endothelial, and selected stromal compartments, where it regulates cell adhesion, polarity, migration, proliferation, barrier organization, and tissue architecture.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com


This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87904-1-RR (EPHA2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.

Immunofluorescent analysis of (-20°C Ethanol) fixed NIH/3T3 cells using EPHA2 antibody (87904-1-RR, Clone: 253572H2) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).



Biolayer interferometry (BLI) kinetic assays of 87904-1-RR against Mouse EPHA2 were performed. The affinity constant is below 1 pM.