

For Research Use Only

STAC3 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87924-1-PBS



Basic Information

Catalog Number:

87924-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5527

GenBank Accession Number:

NM_145064.3

GeneID (NCBI):

246329

UNIPROT ID:

Q96MF2-1

Full Name:

SH3 and cysteine rich domain 3

Calculated MW:

42 kDa

Observed MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

260138A8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87924-1-RR (STAC3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87924-1-PBS in a different storage buffer formulation.