

For Research Use Only

STAC3 Recombinant monoclonal antibody

Catalog Number: 87924-1-RR



Basic Information

Catalog Number: 87924-1-RR	GenBank Accession Number: NM_145064.3	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 246329	CloneNo.: 260138A8
Source: Rabbit	UNIPROT ID: Q96MF2-1	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: SH3 and cysteine rich domain 3	
Immunogen Catalog Number: EG5527	Calculated MW: 42 kDa	
	Observed MW: 55 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls:
Species Specificity: human, mouse, pig	WB : human skeletal muscle tissue, mouse skeletal muscle tissue, pig skeletal muscle tissue

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87924-1-RR (STAC3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.