

For Research Use Only

CXCL13/BCA1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87928-3-PBS



Basic Information

Catalog Number: 87928-3-PBS	GenBank Accession Number: NM_006419.2	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10563	CloneNo.: 253456A10
Source: Rabbit	UNIPROT ID: Q53X90	
Isotype: IgG	Full Name: chemokine (C-X-C motif) ligand 13	
Immunogen Catalog Number: EG7258	Calculated MW: 13 kDa	
	Observed MW: 13 kDa	

Applications

Tested Applications:
WB, IHC, Indirect ELISA

Species Specificity:
human

Background Information

CXCL13 (C-X-C motif chemokine ligand 13), also known as B cell-attracting chemokine 1 (BCA-1) or B-lymphocyte chemoattractant (BLC), is a small cytokine belonging to the CXC chemokine family. It was independently discovered by two research groups in 1998-Legler et al. identified it as BCA-1, while Gunn et al. described it as BLC. CXCL13 is primarily secreted by follicular dendritic cells (FDCs) within B-cell follicles of secondary lymphoid organs, including lymph nodes, spleen, and Peyer's patches.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

