

CXCL13/BCA1 Recombinant monoclonal antibody

Catalog Number: 87928-3-RR

Basic Information

Catalog Number: 87928-3-RR	GenBank Accession Number: NM_006419.2	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 10563	CloneNo.: 253456A10
Source: Rabbit	UNIPROT ID: Q53X90	Recommended Dilutions: WB: 1:500-1:2000 IHC: 1:250-1:1000
Isotype: IgG	Full Name: chemokine (C-X-C motif) ligand 13	
Immunogen Catalog Number: EG7258	Calculated MW: 13 kDa	
	Observed MW: 13 kDa	

Applications

Tested Applications: WB, IHC, ELISA	Positive Controls:
Species Specificity: human	WB : TT cells,
	IHC : human tonsillitis tissue,

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

CXCL13 (C-X-C motif chemokine ligand 13), also known as B cell-attracting chemokine 1 (BCA-1) or B-lymphocyte chemoattractant (BLC), is a small cytokine belonging to the CXC chemokine family. It was independently discovered by two research groups in 1998-Legler et al. identified it as BCA-1, while Gunn et al. described it as BLC. CXCL13 is primarily secreted by follicular dendritic cells (FDCs) within B-cell follicles of secondary lymphoid organs, including lymph nodes, spleen, and Peyer's patches.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



TT cells were subjected to SDS PAGE followed by western blot with 87928-3-RR (CXCL13 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 87928-3-RR (CXCL13 antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 87928-3-RR (CXCL13 antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).