

For Research Use Only

MGA Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87932-2-PBS



Basic Information

Catalog Number: 87932-2-PBS	GenBank Accession Number: NM_001365693.1	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 8972	CloneNo.: 253246G10
Source: Rabbit	UNIPROT ID: O43451-2	
Isotype: IgG	Full Name: maltase-glucoamylase (alpha-glucosidase)	
Immunogen Catalog Number: EG6430	Calculated MW: 312 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

87932-2-PBS targets MGA as part of a matched antibody pair:

MP03262-1: 87932-2-PBS capture and 87932-1-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MGAM concentration was determined to be 439.00 ng/mL in HL-60 cell extract based on a 1.5 mg/mL extract load, 461.37 ng/mL in COLO 320 cell extract based on a 1.2 mg/mL extract load, 144.04 ng/mL in A2780 cell extract based on a 1.2 mg/mL extract load, 173.49 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load and 118.72 ng/mL in Ramos cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP03262-1, Human MGAM Recombinant Matched Antibody Pair - PBS only. 87932-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg6430. 87932-1-PBS was HRP conjugated as the detection antibody. Range: 0.625-40 ng/mL