

For Research Use Only

TFAP2A/AP-2 Recombinant monoclonal antibody

Catalog Number: 87936-1-RR



Basic Information

Catalog Number: 87936-1-RR	GenBank Accession Number: BC017754	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 7020	CloneNo.: 253521E1
Source: Rabbit	UNIPROT ID: P05549	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: transcription factor AP-2 alpha (activating enhancer binding protein 2 alpha)	
Immunogen Catalog Number: AG4112	Calculated MW: 431 aa, 47 kDa	
	Observed MW: 47 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells, BT-549 cells, HeLa cells
Species Specificity: human	

Background Information

The activator protein-2 (AP-2) family of transcription factors comprises five 52-kDa isoforms (AP-2α, AP-2β, AP-2γ, AP-2δ, and AP-2ε), which share a common structure: a proline/glutamine-rich transactivation domain in the N-terminal region and a helix-span-helix domain in the C-terminal region, which mediates dimerization and site-specific DNA binding. Depending on the cellular context, the AP-2 transcription factors are individually associated either with cell differentiation and development or with cancer progression/regression (PMID: 21966377). TFAP2A (AP-2-alpha) is the only AP-2 protein required for early morphogenesis of the lens vesicle. Together with the CITED2 coactivator, stimulates the PITX2 P1 promoter transcription activation (PMID: 11694877).

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87936-1-RR (TFAP2A,AP-2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.