

For Research Use Only

# TCEAL7 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87945-1-PBS



## Basic Information

<b>Catalog Number:</b> 87945-1-PBS	<b>GenBank Accession Number:</b> BC016786	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 56849	<b>CloneNo.:</b> 253550E9
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9BRU2	
<b>Isotype:</b> IgG	<b>Full Name:</b> transcription elongation factor A (SII)-like 7	
<b>Immunogen Catalog Number:</b> EG5177	<b>Calculated MW:</b> 12 kDa	
	<b>Observed MW:</b> 12 kDa	

## Applications

**Tested Applications:**  
WB, Indirect ELISA

**Species Specificity:**  
human

## Background Information

TCEAL7 is a transcription regulatory factor that has a role in the negative regulation of NF-kappa-B signaling at the basal level by regulating transcriptional activity of NF-kappa-B on its target gene promoters. It can associate with cyclin D1 promoter containing Myc E-box sequence and transcriptionally represses cyclin D1 expression. In addition, TCEAL7 regulates telomerase reverse transcriptase expression and telomerase activity in both ALT (alternative lengthening of telomeres) and telomerase-positive cell lines. This antibody can detect a band around 12 kDa, which corresponds to the calculated molecular weight.

## Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87945-1-RR (TCEAL7 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87945-1-PBS in a different storage buffer formulation.