

For Research Use Only

PSG5 Recombinant monoclonal antibody, PBS Only

Catalog Number:87956-1-PBS



Basic Information

Catalog Number:

87956-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG8786

GenBank Accession Number:

NM_002781.4

GeneID (NCBI):

5673

UNIPROT ID:

Q15238

Full Name:

pregnancy specific beta-1-glycoprotein 5

Calculated MW:

38 kDa

Observed MW:

35-45 kDa

Purification Method:

Protein A purification

CloneNo.:

260012B6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87956-1-RR (PSG5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87956-1-PBS in a different storage buffer formulation.