

For Research Use Only

PSG5 Recombinant monoclonal antibody

Catalog Number:87956-1-RR



Basic Information

Catalog Number:

87956-1-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG8786

GenBank Accession Number:

NM_002781.4

GeneID (NCBI):

5673

UNIPROT ID:

Q15238

Full Name:

pregnancy specific beta-1-glycoprotein 5

Calculated MW:

38 kDa

Observed MW:

35-45 kDa

Purification Method:

Protein A purification

CloneNo.:

260012B6

Recommended Dilutions:

WB: 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human placenta tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87956-1-RR (PSG5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.