

For Research Use Only

Mouse FGRL1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88000-1-PBS



Basic Information

Catalog Number: 88000-1-PBS	GenBank Accession Number: NM_001425814.1	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 116701	CloneNo.: 260184E1
Source: Rabbit	UNIPROT ID: Q91V87-1	
Isotype: IgG	Full Name: fibroblast growth factor receptor-like 1	
Immunogen Catalog Number: EG9200	Calculated MW: 57 kDa	
	Observed MW: 50 kDa, 60 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
mouse, rat

Background Information

FGR-like 1 (FGRL1) is a poorly characterized protein-coding gene in mice, annotated based on sequence similarity to members of the Src family kinase-related proteins, particularly FGR, although it is not confirmed to possess canonical kinase activity. Current data suggest that Fgrl1 expression is low and tissue-restricted, with detectable transcripts reported in immune-related tissues and select epithelial compartments, but without a clearly defined cellular distribution.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88000-1-RR (FGRL1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88000-1-PBS in a different storage buffer formulation.