

For Research Use Only

CFAP45 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:88026-2-PBS



Basic Information

Catalog Number: 88026-2-PBS	GenBank Accession Number: BC027874	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 25790	CloneNo.: 260397B9
Source: Rabbit	UNIPROT ID: Q9UL16	
Isotype: IgG	Full Name: coiled-coil domain containing 19	
Immunogen Catalog Number: AG16633	Calculated MW: 551 aa, 66 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

88026-2-PBS targets CFAP45 as part of a matched antibody pair:

MP03337-1: 88026-3-PBS capture and 88026-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CFAP45 concentration was determined to be 85.2 ng/mL in mouse testis tissue extract based on a 4.2 mg/mL extract load, mouse spleen tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03337-1, CFAP45 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88026-3-PBS. Detection antibody: 88026-2-PBS. Standard: Ag16633. Range: 0.078-10 ng/mL.