

For Research Use Only

PSG9 Recombinant monoclonal antibody

Catalog Number: 88037-1-RR



Basic Information

Catalog Number: 88037-1-RR	GenBank Accession Number: NM_002784.5	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 5678	CloneNo.: 260068A8
Source: Rabbit	UNIPROT ID: Q00887-1	Recommended Dilutions: WB: 1:1000-1:4000
Isotype: IgG	Full Name: pregnancy specific beta-1-glycoprotein 9	
Immunogen Catalog Number: EG8782	Calculated MW: 48 kDa	
	Observed MW: 45-70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human placenta tissue,
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88037-1-RR (PSG9 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.