

For Research Use Only

RNASE2 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88054-1-PBS



Basic Information

Catalog Number:

88054-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28325

GenBank Accession Number:

BC093678

GeneID (NCBI):

6036

UNIPROT ID:

P10153

Full Name:

ribonuclease, RNase A family, 2 (liver, eosinophil-derived neurotoxin)

Calculated MW:

161 aa, 18 kDa

Observed MW:

18 kDa

Purification Method:

Protein A purification

CloneNo.:

260304C6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human urine sample were subjected to SDS PAGE followed by western blot with 88054-1-RR (RNASE2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88054-1-PBS in a different storage buffer formulation.