

For Research Use Only

CSTF2T Recombinant monoclonal antibody, PBS Only

Catalog Number: 88063-1-PBS



Basic Information

Catalog Number:

88063-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5411

GenBank Accession Number:

BC028239

GeneID (NCBI):

23283

UNIPROT ID:

Q9HOL4

Full Name:

cleavage stimulation factor, 3' pre-RNA, subunit 2, 64kDa, tau variant

Calculated MW:

64 kDa

Observed MW:

~70 kDa

Purification Method:

Protein A purification

CloneNo.:

260142B6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88063-1-RR (CSTF2T antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88063-1-PBS in a different storage buffer formulation.



rat brain tissue were subjected to SDS PAGE followed by western blot with 88063-1-RR (CSTF2T antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88063-1-PBS in a different storage buffer formulation.