

For Research Use Only

Integrin alpha-9 Recombinant monoclonal antibody

Catalog Number: 88081-2-RR



Basic Information

Catalog Number: 88081-2-RR	GenBank Accession Number: NM_002207.2	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 3680	CloneNo.: 260170E4
Source: Rabbit	UNIPROT ID: Q13797	Recommended Dilutions: WB: 1:2000-1:10000
Isotype: IgG	Full Name: integrin, alpha 9	
Immunogen Catalog Number: EG6174	Calculated MW: 114 kDa	
	Observed MW: 140-150 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A-673 cells, A549 cells, HeLa cells, HepG2 cells
Species Specificity: human	

Background Information

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated α and β subunits (PMID: 9779984). Integrin alpha 9 contains a large N-terminal extracellular domain with 7 conserved repeats of putative metal-binding domains, a transmembrane segment, and a short C-terminal cytoplasmic tail (PMID: 8245132). Integrin alpha 9 interacts with integrin beta 1, generating $\alpha 9\beta 1$ heterodimer, which is expressed in many cell types, including epithelial cells, neutrophils, hepatocytes, muscle, and endothelial cells (PMID: 21975548). It binds to a diversity of ligands in the extracellular matrix, like the α -disintegrin and metalloprotease (ADAM) family, THBS1, EMILIN1, VEGF, tenascin-C, osteopontin, fibronectin, and VCAM1 (PMID: 21975548; 33790909). Integrin $\alpha 9\beta 1$ has been implicated in diverse biological functions, including angiogenesis, lymphangiogenesis, lymphatic valve morphogenesis, granulopoiesis and tumorigenesis (PMID: 20179422).

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88081-2-RR (Integrin alpha-9 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.