

For Research Use Only

PODO Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88093-2-PBS



Basic Information

Catalog Number: 88093-2-PBS	GenBank Accession Number: NM_014625.4	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7827	CloneNo.: 260462A5
Source: Rabbit	UNIPROT ID: Q9NP85-1	
Isotype: IgG	Full Name: nephrosis 2, idiopathic, steroid-resistant (podocin)	
Immunogen Catalog Number: EG9317	Calculated MW: 42 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

88093-2-PBS targets PODO as part of a matched antibody pair:

MP03429-1: 88093-3-PBS capture and 88093-2-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean NPHS2 concentration was determined to be 7.6 ng/mL in mouse kidney extract based on a 8.4 mg/mL extract load and 0.5 ng/mL in mouse brain extract based on a 6.3 mg/mL extract load.



Sandwich ELISA standard curve of MP03429-1, Human NPHS2 Recombinant Matched Antibody Pair - PBS only. 88093-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg9317. 88093-2-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL