

Nectin-1/PVRL1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88101-1-PBS

Basic Information

Catalog Number: 88101-1-PBS	GenBank Accession Number: NM_002855.5	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5818	CloneNo.: 260037D11
Source: Rabbit	UNIPROT ID: Q15223-1	
Isotype: IgG	Full Name: poliovirus receptor-related 1 (herpesvirus entry mediator C)	
Immunogen Catalog Number: EG1707	Calculated MW: 57 kDa	
	Observed MW: 90 kDa, 57 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, rat

Background Information

Nectins are a family of immunoglobulin-like cell adhesion molecules comprising four members, nectin-1 to nectin-4, each consisting of an extracellular region with three immunoglobulin-like loops, a single transmembrane domain and a cytoplasmic domain (PMID: 12901789, 18703497). Nectins are involved in cell adhesion, migration, polarization, morphogenesis, proliferation, differentiation, and survival. Nectin-1, also known as PVRL1, HVEC and CD111, is widely expressed in a variety of cell types including epithelial cells, neurons, and fibroblasts (PMID: 12901789). It promotes cell-cell contacts by forming homophilic or heterophilic trans-dimers. Nectin-1 acts as a receptor for several alphaherpesviruses, including herpes simplex viruses (HSV) 1 and 2, porcine pseudorabies virus (PRV), and bovine herpesvirus 1 (BHV-1) (PMID: 9616127). Nectin-1 is a glycoprotein, the glycosylated form has a larger apparent molecular weight than the calculated molecular weight of 57 kDa (PMID: 12376527; 15731277).

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88101-1-RR (Nectin-1/PVRL1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88101-1-PBS in a different storage buffer formulation.



rat cerebellum tissue were subjected to SDS PAGE followed by western blot with 88101-1-RR (Nectin-1/PVRL1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88101-1-PBS in a different storage buffer formulation.