

For Research Use Only

PEBP1/RKIP Recombinant monoclonal antibody, PBS Only

Catalog Number: 88106-1-PBS



Basic Information

Catalog Number: 88106-1-PBS	GenBank Accession Number: BC008714	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5037	CloneNo.: 260328F7
Source: Rabbit	UNIPROT ID: P30086	
Isotype: IgG	Full Name: phosphatidylethanolamine binding protein 1	
Immunogen Catalog Number: AG0864	Calculated MW: 21 kDa	
	Observed MW: 23 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse

Background Information

Raf kinase inhibitor protein (RKIP), also termed phosphatidylethanolamine binding protein PEBP1, was initially identified as a potent inhibitor of Raf-1/MEK/ERK, NF- κ B, and G-protein-coupled receptor signaling pathways. Later RKIP has been further identified as a metastasis suppressor and its loss of expression has been reported in various cancers. Its expression has been proposed as a prognostic marker for patients diagnosed with the above cancers.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88106-1-RR (PEBP1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88106-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 88106-1-RR (PEBP1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88106-1-PBS in a different storage buffer formulation.