

For Research Use Only

SH3RF1/POSH Recombinant monoclonal antibody, PBS Only

Catalog Number: 88135-1-PBS



Basic Information

Catalog Number:

88135-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6243

GenBank Accession Number:

BC053671

GeneID (NCBI):

57630

UNIPROT ID:

Q7Z6J0

Full Name:

SH3 domain containing ring finger 1

Calculated MW:

93 kDa

Observed MW:

93-100 kDa

Purification Method:

Protein A purification

CloneNo.:

260354A12

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

POSH (plenty of SH3 domains) is a multidomain GTP-Rac1-interacting protein that acts as a scaffold for JNK activation, leading to apoptosis in neuronal and non-neuronal cells. It consists a RING-finger domain that possesses ubiquitin-ligase activity and four SH3 domains.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88135-1-RR (SH3RF1/POSH antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88135-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 88135-1-RR (SH3RF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88135-1-PBS in a different storage buffer formulation.