

For Research Use Only

SH3RF1/POSH Recombinant monoclonal antibody

Catalog Number: 88135-1-RR



Basic Information

Catalog Number: 88135-1-RR	GenBank Accession Number: BC053671	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 57630	CloneNo.: 260354A12
Source: Rabbit	UNIPROT ID: Q7Z6J0	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: SH3 domain containing ring finger 1	
Immunogen Catalog Number: AG6243	Calculated MW: 93 kDa	
	Observed MW: 93-100 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HEK-293 cells, DU 145 cells, A431 cells, Neuro-2a cells, C6 cells, mouse brain tissue
Species Specificity: human, mouse, rat	

Background Information

POSH (plenty of SH3 domains) is a multidomain GTP-Rac1-interacting protein that acts as a scaffold for JNK activation, leading to apoptosis in neuronal and non-neuronal cells. It consists a RING-finger domain that possesses ubiquitin-ligase activity and four SH3 domains.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88135-1-RR (SH3RF1/POSH antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 88135-1-RR (SH3RF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.