

For Research Use Only

ROBO4 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88139-2-PBS



Basic Information

Catalog Number: 88139-2-PBS	GenBank Accession Number: NML_019055.6	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 54538	CloneNo.: 260378E9
Source: Rabbit	UNIPROT ID: Q8WZ75-1	
Isotype: IgG	Full Name: roundabout homolog 4, magic roundabout (Drosophila)	
Immunogen Catalog Number: EG7821	Calculated MW: 107 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

88139-2-PBS targets ROBO4 as part of a matched antibody pair:

MP03363-1: 88139-3-PBS capture and 88139-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ROBO4 concentration was determined to be 94.5 ng/mL in Human plasma 1, 74.4 ng/mL in Human plasma 2.



Cytometric bead array standard curve of MP03363-1, ROBO4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88139-3-PBS. Detection antibody: 88139-2-PBS. Standard: Eg7821. Range: 0.312-40 ng/mL.