

DHX9 Recombinant monoclonal antibody

Catalog Number: 88165-1-RR

Basic Information

Catalog Number: 88165-1-RR	GenBank Accession Number: BC014246	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 1660	CloneNo.: 260493F8
Source: Rabbit	UNIPROT ID: Q08211	Recommended Dilutions: WB: 1:5000-1:50000 IHC: 1:500-1:2000
Isotype: IgG	Full Name: DEAH (Asp-Glu-Ala-His) box polypeptide 9	
Immunogen Catalog Number: AG12807	Calculated MW: 1270 aa, 141 kDa Observed MW: 140 kDa	

Applications

Tested Applications:

WB, IHC, ELISA

Species Specificity:

human, mouse, rat

Note-IHC: suggested antigen retrieval with *TE buffer pH 9.0*; (*) Alternatively, antigen retrieval may be performed with *citrate buffer pH 6.0*

Positive Controls:

WB : HeLa cells, HepG2 cells, THP-1 cells, K-562 cells, HEK-293 cells, mouse testis tissue

IHC : human skin tissue,

Background Information

RNA helicases play important roles in transcription, RNA processing, translation, and RNA replication. DEAD box proteins are putative RNA helicases that have a characteristic Asp-Glu-Ala-Asp (DEAD) box as 1 of 8 highly conserved sequence motifs. DHX9 a member of the DEAH family of proteins, which possess a double-stranded RNA-binding domain (dsRBD) and a helicase domain (PMID:20569003). It unwinds double-stranded DNA and RNA in a 3' to 5' direction. Alteration of secondary structure of DHX9 may subsequently influence interactions with proteins or other nucleic acids. It is also a component of the CRD-mediated complex that promotes MYC mRNA stability. In addition, it is involved with LARP6 in the stabilization of type I collagen mRNAs for CO1A1 and CO1A2 (PMID: 19029303, 22190748).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88165-1-RR (DHX9 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human skin tissue slide using 88165-1-RR (DHX9 antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human skin tissue slide using 88165-1-RR (DHX9 antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).