

For Research Use Only

Cleaved GSDME (Asp270), C-terminal Recombinant monoclonal antibody, PBS Only



Catalog Number: 88249-1-PBS

Basic Information

Catalog Number: 88249-1-PBS	GenBank Accession Number: BC019689	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 1687	CloneNo.: 253373H3
Source: Rabbit	UNIPROT ID: O60443	
Isotype: IgG	Full Name: deafness, autosomal dominant 5	
	Calculated MW: 496 aa, 55 kDa	
	Observed MW: 25 kDa	

Applications

Tested Applications:
WB, FC (Intra), Indirect ELISA

Species Specificity:
human

Background Information

GSDME, also known as DFNA5 (deafness, autosomal dominant 5), is a member of the gasdermin family of pore-forming proteins. In its full-length form, GSDME is autoinhibited. Upon apoptotic stimulation, active caspase-3 specifically cleaves GSDME at a conserved aspartate residue, Asp270 (D270). This cleavage of Gasdermin E separates the N-terminal pore-forming domain from the C-terminal repressor domain. The released N-terminal fragment then inserts into the plasma membrane, forms large pores, and induces a lytic, inflammatory form of cell death known as pyroptosis. Therefore, Asp270 cleavage is a critical molecular switch that converts caspase-3-mediated apoptosis into pyroptosis, linking GSDME/DFNA5 to cancer therapy, chemotherapy-induced toxicity, and hearing loss pathology.

Storage



Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Doxorubicin hydrochloride treated A431 cells were subjected to SDS PAGE followed by western blot with 88249-1-RR (Cleaved GSDME (Asp270), C-terminal antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88249-1-PBS in a different storage buffer formulation.

1×10^6 SH-SY5Y cells untreated (dashed lines) or treated with etoposide which intracellularly stained with 0.25 ug cleaved GSDME (Asp270), C-terminal Recombinant monoclonal antibody (88249-1-RR, Clone:253373H3) and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with Intracellular Flow Cytometry Fixation & Permeabilization Buffer Kit (PF00019).

