

For Research Use Only

CPA1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88295-1-PBS



Basic Information

Catalog Number: 88295-1-PBS	GenBank Accession Number: BC005279	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 1357	CloneNo.: 260663D6
Source: Rabbit	UNIPROT ID: P15085	
Isotype: IgG	Full Name: carboxypeptidase A1 (pancreatic)	
Immunogen Catalog Number: AG8586	Calculated MW: 419 aa, 47 kDa	
	Observed MW: 44-47 kDa	

Applications

Tested Applications:
WB, IF-P, Indirect ELISA

Species Specificity:
human, mouse, rat, pig

Background Information

CPA1(Carboxypeptidase A1) is also named as CPA and belongs to the the peptidase M14 family. CPA1 is one of the genes whose expression in the pancreas and the protein which it encodes is synthesized as an inactive precursor, proCPA1, which is processed to the active enzyme by the proteolytic removal of the 95-amino acid N-terminal prodomain. The mature CPA1 is not secreted and is trapped and degraded intracellularly(PMID:8639628). This antibody is speicific to CPA1.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.



Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 88295-1-RR (CPA1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88295-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 88295-1-RR (CPA1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88295-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded rat pancreas tissue using CPA1 antibody (88295-1-RR, Clone: 260663D6) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 88295-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse pancreas tissue using CPA1 antibody (88295-1-RR, Clone: 260663D6) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 88295-1-PBS in a different storage buffer formulation.