

For Research Use Only

Sam68 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88301-3-PBS



Basic Information

Catalog Number: 88301-3-PBS	GenBank Accession Number: BC000717	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10657	CloneNo.: 260478D5
Source: Rabbit	UNIPROT ID: Q07666	
Isotype: IgG	Full Name: KH domain containing, RNA binding, signal transduction associated 1	
Immunogen Catalog Number: AG0235	Calculated MW: 48 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

88301-3-PBS targets Sam68 as part of a matched antibody pair:

MP03453-2: 88301-1-PBS capture and 88301-3-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Sam68 concentration was determined to be 17.46 ng/mL in A2780 cell extract based on a 1.2 mg/mL extract load, 6.52 ng/mL in Caco-2 cell extract based on a 1.2 mg/mL extract load and 683.88 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP03453-2, Human Sam68 Recombinant Matched Antibody Pair - PBS only. 88301-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0235. 88301-3-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL