

For Research Use Only

CARNS1 Recombinant monoclonal antibody

Catalog Number: 88342-1-RR



Basic Information

Catalog Number: 88342-1-RR	GenBank Accession Number: BC036557	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 57571	CloneNo.: 260860B7
Source: Rabbit	UNIPROT ID: A5YM72	Recommended Dilutions: WB: 1:2000-1:10000
Isotype: IgG	Full Name: ATP-grasp domain containing 1	
Immunogen Catalog Number: AG22761	Calculated MW: 827 aa, 88 kDa	
	Observed MW: 88-103 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : pig skeletal muscle tissue, human skeletal muscle tissue
Species Specificity: human	

Background Information

CARNS1, located at human chromosome 11q13.2, encodes carnosine synthase 1 belonging to the ATP-grasp ATPase family (EC 6.3.2.11). It mediates ATP-dependent synthesis of carnosine in skeletal muscle and homocarnosine in the central nervous system, two vital histidine-containing dipeptides with antioxidant and pH-buffering effects. This gene is highly conserved across vertebrates and preferentially expressed in brain and muscle tissues. Aberrant CARNS1 expression disturbs dipeptide homeostasis, linking to aging-related sensory decline and tumor progression via impaired carnosine-mediated tumor suppression. (PMID: 20097752)

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88342-1-RR (ATPGD1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.