

For Research Use Only

EIF4B Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88367-2-PBS



Basic Information

Catalog Number: 88367-2-PBS	GenBank Accession Number: BC073154	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 1975	CloneNo.: 260892B9
Source: Rabbit	UNIPROT ID: P23588	
Isotype: IgG	Full Name: eukaryotic translation initiation factor 4B	
Immunogen Catalog Number: AG12325	Calculated MW: 611 aa, 69 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

88367-2-PBS targets EIF4B as part of a matched antibody pair:

MP03486-1: 88367-3-PBS capture and 88367-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean EIF4B concentration was determined to be 247.9 ng/mL in HeLa cell extract based on a 1.0 mg/mL extract load, 122.5 ng/mL in HEK-293 cell extract based on a 1.0 mg/mL extract load, 87.1 ng/mL in MCF-7 cell extract based on a 1.0 mg/mL extract load.



Cytometric bead array standard curve of MP03486-1, EIF4B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88367-3-PBS. Detection antibody: 88367-2-PBS. Standard: Ag12325. Range: 0.156-20 ng/mL.