

For Research Use Only

ZCCHC6 Recombinant monoclonal antibody

Catalog Number: 88391-1-RR



Basic Information

Catalog Number: 88391-1-RR	GenBank Accession Number: BC032456	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 79670	CloneNo.: 260850G10
Source: Rabbit	UNIPROT ID: Q5VYS8	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: zinc finger, CCHC domain containing 6	
Immunogen Catalog Number: AG18240	Calculated MW: 1495 aa, 171 kDa	
	Observed MW: 171 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HEK-293T cells, MCF-7 cells, PC-3 cells, mouse brain tissue
Species Specificity: human, mouse	

Background Information

ZCCHC6 also named as KIAA1711 ,HS2 or TUT7 is a 1495 amino acid protein, which contains three CCHC-type zinc fingers and two PAP-associated domains. A member of the DNA polymerase type-B-like family, ZCCHC6 as a uridylyltransferase that mediates RNA uridylation.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88391-1-RR (ZCCHC6 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.